

# Anti-Pdgfrb antibody



Catalog Number: 176523

## Product name

Anti-Pdgfrb antibody

## Specificity

Human, Rat

## Antibody description

Rabbit monoclonal antibody to Pdgfrb

## Preparation

This antigen of this antibody was synthetic phospho-peptide corresponding to residues surrounding tyr740 of human pdgfr beta.

## Formulation

Liquid, 1\*TBS (pH7.4), 0.05% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications

WB, ICC/IF

## Dilutions

WB: 1:1,000

ICC/IF: 1:50-1:200

## Validations




Fig1: ICC staining PDGF Receptor beta(phospho Y740) in C2C12 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig2: ICC staining PDGF Receptor beta(phospho Y740) in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig3: ICC staining PDGF Receptor beta(phospho Y740) in L6 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in

# Anti-Pdgfrb antibody



Catalog Number: 176523

---

paraformaldehyde, permeabilised with 0.25%

Triton X100/PBS.