

# Anti-MROH7 antibody

Catalog Number: 176651

## General Information

### Product name:

Anti-MROH7 antibody

### Immunogen:

### Applications:

WB, IHC-P

### Specificity:

Human

### Antibody description:

Mouse monoclonal antibody to MROH7

## Properties

### Preparation:

This antigen of this antibody was recombinant protein

### Formulation:

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

### Storage:

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

### Clonality:

Monoclonal

**Ig Type:**IgG1

## Applications

WB: 1:500-1:1,000

IHC-P: 1:50-1:200

## Validations




Fig1: Western blot analysis of C1orf175 on recombinant protein using anti-C1orf175 antibody at 1/1,000 dilution.




Fig2: Immunohistochemical analysis of paraffin-embedded human kidney tissue using anti-C1orf175 antibody. Counter stained with hematoxylin.