

# Anti-Lrp1 antibody



Catalog Number: 176655

## General Information

### Product name:

Anti-Lrp1 antibody

### Immunogen:

### Applications:

WB, IHC-P

### Specificity:

Mouse

### Antibody description:

Mouse monoclonal antibody to Lrp1

## Properties

### Preparation:

This antigen of this antibody was recombinant protein

### Formulation:

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

### Storage:

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

### Clonality:

Monoclonal

**Ig Type:** IgG2a

## Applications

WB: 1:1,000

IHC-P: 1:50

## Validations




Fig1: Western blot analysis on recombinant protein using anti-LRP-1 Mouse mAb (Cat. # 176655#).




Fig2: Immunohistochemical analysis of paraffin-embedded mouse liver tissue using anti-LRP-1 Mouse mAb (Cat. # 176655#).