

Anti-IL31 antibody



Catalog Number: 176656

General Information

Product name:

Anti-IL31 antibody

Immunogen:

Applications:

WB, ELISA

Specificity:

Human

Antibody description:

Mouse monoclonal antibody to IL31

Properties

Preparation:

This antigen of this antibody was recombinant protein

Formulation:

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage:

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality:

Monoclonal

Ig Type: IgG2a

Applications

WB: 1:2,000

ELISA: 1:5,000-1:10,000

Validations




Fig1: Western blot analysis on bioactive recombinant human IL-31 protein using anti-IL-31 Mouse mAb (Cat. # 176656#).