

Anti-NXPE3 antibody



Catalog Number: 176658

General Information

Product name:

Anti-NXPE3 antibody

Immunogen:

Applications:

WB, ICC

Specificity:

Human

Antibody description:

Mouse monoclonal antibody to NXPE3

Properties

Preparation:

This antigen of this antibody was recombinant protein

Formulation:

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage:

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality:

Monoclonal

Ig Type: IgG1

Applications

WB: 1:1,000-1:2,000

ICC: 1:50-1:200

Validations




Fig1: Western blot analysis of Fam55C on Fam55C-GST recombinant protein lysates using anti-Fam55C antibody at 1/5,000 dilution.




Fig2: ICC staining Fam55C in NCCIT cells (green). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.