

Anti-nusA antibody



Catalog Number: 176659

General Information

Product name:

Anti-nusA antibody

Immunogen:

Applications:

WB, ELISA

Specificity:

Escherichia coli

Antibody description:

Mouse monoclonal antibody to nusA

Properties

Preparation:

This antigen of this antibody was recombinant protein

Formulation:

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.

Preservative: 0.05% Sodium Azide.

Storage:

Store at -20°C

Clonality:

Monoclonal

Ig Type: IgG1

Applications

WB: 1:2,000-1:5,000

Validations




Fig1: Western blot analysis on E.coli cell lysates using anti- NusA Mouse mAb (Cat. # 176659#).