

Anti-Sox10 antibody

Catalog Number: 176660

General Information

Product name:

Anti-Sox10 antibody

Immunogen:

Applications:

WB, IHC-P

Specificity:

Mouse

Antibody description:

Mouse monoclonal antibody to Sox10

Properties

Preparation:

This antigen of this antibody was recombinant protein

Formulation:

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage:

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality:

Monoclonal

Ig Type: IgG1

Applications

WB: 1:500-1:1000

IHC-P: 1:50

Validations




Fig1: Western blot analysis on F9 (1) and mouse fetal brain (2) cell lysates using anti-SOX10 Mouse mAb.




Fig2: Immunohistochemical analysis of paraffin-embedded mouse fetal brain (13.5 days) tissue using anti- SOX10 Mouse mAb.