

# Anti-UTF1 antibody



Catalog Number: 176665

## General Information

### Product name:

Anti-UTF1 antibody

### Immunogen:

### Applications:

WB, ICC

### Specificity:

Human, Mouse, Rat

### Antibody description:

Mouse monoclonal antibody to UTF1

## Properties

### Preparation:

This antigen of this antibody was peptide

### Formulation:

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

### Storage:

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

### Clonality:

Monoclonal

### Ig Type: IgG2b

## Applications

WB: 1:500-1:1,000

ICC: 1:500-1:1,000

## Validations




Fig1: Western blot analysis on D3 cell lysates using anti- UTF1 mouse mAb.