

Anti-IGHA1 antibody



Catalog Number: 176673

General Information

Product name:

Anti-IGHA1 antibody

Immunogen:

Applications:

WB, ELISA

Specificity:

Human

Antibody description:

Mouse monoclonal antibody to IGHA1

Properties

Preparation:

This antigen of this antibody was native protein

Formulation:

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage:

Store at +4°C after thawing. Aliquot store at -20°C

or -80°C. Avoid repeated freeze / thaw cycles.

Clonality:

Monoclonal

Ig Type: IgG1

Applications

WB: 1:500-1:2000

ELISA: Use at an assay dependent concentration.

Validations




Fig1:; Western blot analysis of IgA on human plasma tissue lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody (1/500) was used in 5% BSA at room temperature for 2 hours. Goat Anti-Mouse IgG - HRP Secondary Antibody (HA1006) at 1:5,000 dilution was used for 1 hour at room temperature.