

Anti-IGHA1 antibody



Catalog Number: 176675

General Information

Product name:

Anti-IGHA1 antibody

Immunogen:

Applications:

Elisa, WB

Specificity:

Human

Antibody description:

Mouse monoclonal antibody to IGHA1

Properties

Preparation:

This antigen of this antibody was native protein.

Formulation:

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage:

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality:

Monoclonal

Ig Type: IgG2a

Applications

Elisa: 1:5000

WB: 1:500-1:2,000

Validations




Fig1: Western blot analysis of IgA on human plasma lysate. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody was used at a 1:500 dilution in 5% BSA at room temperature for 2 hours. Goat Anti-Mouse IgG - HRP Secondary Antibody (HA1006) at 1:5,000 dilution was used for 1 hour at room temperature.