

DM1 Monoclonal Antibody



Catalog Number: 176726

Product name

DM1 Monoclonal Antibody

Specificity

Chemical

Preparation

DM1 conjugated with BSA, Protein G affinity purified

Formulation

PBS(pH7.4)

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Mouse IgG2b, κ

Ig Type

Human IgG


Applications

ELISA

Dilutions


ELISA: 1:5,000-1:20,000

Validations



ELISA Binding Assay of Anti-DM1 Antibody (MM56B) to different proteins.

ELISA analysis of Anti-DM1 Monoclonal Antibody with DM1-BSA, DM1-KLH, BSA, KLH and Kadcylla. Anti-DM1 Monoclonal Antibody can bind to DM1-BSA, DM1-KLH or DM1 ADC but not carrier



proteins.

ELISA Binding Assay of Anti-DM1 Antibody (MM56B)

Immobilized DM1 ADC at 1µg/mL (100µL/well) can bind anti-DM1 mouse monoclonal antibody.