

# DM1 Monoclonal Antibody



Catalog Number: 176726

## Product name

DM1 Monoclonal Antibody

## Specificity

Chemical

## Preparation

DM1 conjugated with BSA, Protein G affinity purified

## Formulation

PBS(pH7.4)

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Mouse IgG2b,  $\kappa$

## Ig Type

Human IgG


## Applications

ELISA

## Dilutions


ELISA: 1:5,000-1:20,000

## Validations



ELISA Binding Assay of Anti-DM1 Antibody (MM56B) to different proteins.

ELISA analysis of Anti-DM1 Monoclonal Antibody with DM1-BSA, DM1-KLH, BSA, KLH and Kadcyla. Anti-DM1 Monoclonal Antibody can bind to DM1-BSA, DM1-KLH or DM1 ADC but not carrier



proteins.

ELISA Binding Assay of Anti-DM1 Antibody (MM56B)

Immobilized DM1 ADC at 1  $\mu$ g/mL (100  $\mu$ L/well) can bind anti-DM1 mouse monoclonal antibody.