

# MMAE Monoclonal Antibody



Catalog Number: 176727

## Product name

MMAE Monoclonal Antibody

## Specificity

Chemical

## Antibody description

Product has been tested in direct ELISA.

## Preparation

MMAE conjugated with KLH. Protein G affinity purified

## Formulation

PBS(pH7.4)

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1,  $\kappa$


## Applications

ELISA

## Dilutions


ELISA: 1:5,000-1:20,000

## Validations



ELISA Binding Assay of Anti-MMAE Antibody (MM508H) to different proteins.

ELISA analysis of Anti-MMAE Monoclonal Antibody with MMAE-BSA, MMAE ADC and BSA. Anti-MMAE Monoclonal Antibody can bind to MMAE-BSA or MMAE ADC but not carrier protein.



ELISA Binding Assay of Anti-MMAE Antibody (MM508H)

Immobilized MMAE ADC at 2 $\mu$ g/mL (100 $\mu$ L/well) can bind anti-MMAE mouse monoclonal antibody.