

Anti-MMAF Monoclonal Antibody



Catalog Number: 176728

Product name

ELISA: 1:5,000-1:20,000

Anti-MMAF Monoclonal Antibody

Specificity

Chemical

Antibody description

Product has been tested in direct ELISA.

Preparation

MMAF conjugated with BSA. Protein G affinity purified

Formulation

PBS(pH7.4)

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG2b, κ


Applications

ELISA

Dilutions

Validations


Mouse anti-MMAF Antibody ELISA
0.1 μ g/ml



ELISA Binding Assay of Anti-MMAF Antibody (MM734C) to different proteins.

ELISA analysis of Anti-MMAF Monoclonal Antibody with MMAF-BSA, MMAF-KLH, BSA and KLH. Anti-MMAF Monoclonal Antibody can bind to MMAF-BSA or MMAF-KLH, but not carrier

anti-MMAF mouse mAb MM734C



proteins.

ELISA Binding Assay of Anti-MMAF Antibody (MM734C).

Immobilized MMAF ADC at 2 μ g/mL (100 μ L/well) can bind anti-MMAF mouse monoclonal antibody.