

# Tavolixizumab (Synonyms: TAVOLIMAB [USAN]; Tavolimab [WHO-DD])



Catalog Number: 176802

## Product name

Tavolixizumab (Synonyms: TAVOLIMAB [USAN]; Tavolimab [WHO-DD])

## Specificity

Human TNFRSF4 / OX40 / CD134

## Antibody description

Anti-TNFRSF4 / OX40 / CD134 Antibody (Tavolixizumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-TNFRSF4 / OX40 / CD134 Antibody (tavolixizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-TNFRSF4 / OX40 / CD134 Antibody (tavolixizumab) is 97.98%, determined



by SEC-HPLC.

Bioactivity: ELISA

Immobilized human CD134 His at 2 µg/mL can bind Anti-TNFRSF4 / OX40 / CD134 Antibody (tavolixizumab), EC50=0.03236 µg/mL.