

# Divozilimab reference antibody



Catalog Number: 176808

## Product name

Divozilimab reference antibody

## Specificity

Human CD20 / MS4A1

## Antibody description

Anti-CD20 / MS4A1 Antibody (Divozilimab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-CD20 Antibody (divozilimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-CD20 Antibody (divozilimab) is more than 95%, determined by SEC-HPLC.



## ELISA

Immobilized Human CD20 VLP at 5 ug/mL (30 ul/well) can bind Anti-CD20 / MS4A1 Antibody (divozilimab), The EC50 was approximately 0.01698 ug/ml.