

Milatumzumab reference antibody



Catalog Number: 176812

Product name

Milatumzumab reference antibody

Specificity

Human CD74

Antibody description

Anti-CD74 Antibody (Milatumzumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations




SDS-PAGE

Anti-CD74 Antibody (milatumzumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-CD74 Antibody (milatumzumab) is 99.31%, determined by SEC-HPLC.



ELISA

Immobilized human CD74 Fc at 2 ug/mL can bind Anti-CD74 Antibody (milatumzumab), EC50=0.04337 ug/mL.