

Tremelimumab (Synonyms: TICILIMUMAB)



Catalog Number: 176824

Product name

Tremelimumab (Synonyms: TICILIMUMAB)

Specificity

Human CTLA-4 / CD152

Antibody description

Anti-CTLA-4 / CD152 Antibody (Tremelimumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations




SDS-PAGE

Anti-CTLA-4 / CD152 Antibody (tremelimumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-CTLA-4 / CD152 Antibody (tremelimumab) is 100%, determined by SEC-



HPLC.

Bioactivity: ELISA

Immobilized human CTLA 4 His at 2 ug/mL can bind Anti-CTLA-4 / CD152 Antibody (tremelimumab), EC50=0.004377 ug/mL.