

# Emibetuzumab (Synonyms: LA-480)



Catalog Number: 176827

## Product name

Emibetuzumab (Synonyms: LA-480)

## Specificity

Human HGFR / c-Met

## Antibody description

Anti-HGFR / c-Met Antibody (Emibetuzumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG4SP

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-HGFR / c-Met Antibody (emibetuzumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-HGFR / c-Met Antibody (emibetuzumab) is 99.28%, determined by SEC-



## HPLC.

Bioactivity: ELISA

Immobilized human c MET His at 2 ug/mL can bind Anti-HGFR / c-Met Antibody (emibetuzumab), EC50=0.004797 ug/mL.