

# Andecaliximab reference antibody



Catalog Number: 176828

## Product name

Andecaliximab reference antibody

## Specificity

Human MMP9

## Antibody description

Anti-MMP9 Antibody (Andecaliximab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG4

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-MMP9 Antibody (andecaliximab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



## SEC-HPLC

The purity of Anti-MMP9 Antibody (andecaliximab) is more than 95%, determined by



## SEC-HPLC.

Bioactivity: ELISA

Immobilized human MMP9 His at 2 µg/mL can bind Anti-MMP9 Antibody (andecaliximab), EC50=0.006577 µg/mL.