

Ianalumab reference antibody



Catalog Number: 176835

Product name

Ianalumab reference antibody

Specificity

Human TNFRSF13C / BAFFR / CD268

Antibody description

Anti-TNFRSF13C / BAFFR / CD268 Antibody (Ianalumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations




SDS-PAGE

Anti-TNFRSF13C / BAFFR / CD268 Antibody (ianalumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-TNFRSF13C / BAFFR / CD268 Antibody (ianalumab) is 100%, determined by



SEC-HPLC.

Bioactivity: ELISA

Immobilized human GPA33 / A33, Fc Tag at 2 ug/mL can bind ianalumab, EC50=0.008489 ug/mL.