

# Utomilumab reference antibody



Catalog Number: 176836

## Product name

Utomilumab reference antibody

## Specificity

Human TNFRSF9 / 4-1BB / CD137

## Antibody description

Anti-TNFRSF9 / 4-1BB / CD137 Antibody (Utomilumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG2SA

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-TNFRSF9 / 4-1BB / CD137 Antibody (utomilumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-TNFRSF9 / 4-1BB / CD137 Antibody (utomilumab) is 98.78%, determined by



## SEC-HPLC.

Bioactivity: ELISA

Immobilized human TNFRSF9 / 4-1BB / CD137, His Tag at 2 ug/mL can bind utomilumab, EC50=0.006915 ug/mL.