

# Trastuzumab duocarmazine reference antibody



Catalog Number: 176852

## Product name

Trastuzumab duocarmazine reference antibody

## Specificity

Human ERBB2 / HER2 / CD340

## Antibody description

Anti-ERBB2 / HER2 / CD340 Antibody  
(Trastuzumab duocarmazine)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-ERBB2 / HER2 / CD340 Antibody (Trastuzumab duocarmazine) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-ERBB2 / HER2 / CD340 Antibody (Trastuzumab duocarmazine) is more than 95%, determined by SEC-HPLC.