

# Ruplizumab reference antibody



Catalog Number: 176857

## Product name

Ruplizumab reference antibody

## Specificity

Human TNFSF5 / CD40L / CD154

## Antibody description

Anti-TNFSF5 / CD40L / CD154 Antibody (Ruplizumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-TNFSF5 / CD40L / CD154 Antibody (ruplizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-TNFSF5 / CD40L / CD154 Antibody (ruplizumab) is 99%, determined by SEC-HPLC.