

# Enokizumab (Synonyms: 7F3COM2H2)



Catalog Number: 176884

## Product name

Enokizumab (Synonyms: 7F3COM2H2)

## Specificity

Human IL-9

## Antibody description

Anti-IL-9 Antibody (Enokizumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

1×PBS, pH 6.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-IL-9 Antibody (enokizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



## SEC-HPLC

The purity of Anti-IL-9 Antibody (enokizumab) is more than 95%, determined by SEC-HPLC.