

# Narnatumab reference antibody



Catalog Number: 176885

## Product name

Monoclonal

Narnatumab reference antibody

## Ig Type

## Specificity

IgG1

Human MSPR / RON / CD136

## Applications

## Antibody description

ELISA

Anti-MSPR / RON / CD136 Antibody (Narnatumab)

## Validations

## Preparation

Recombinant expression and purified from CHO cells.


## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality



## SDS-PAGE

Anti-MSPR / RON / CD136 Antibody (narnatumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.