

Osocimab (Synonyms: BAY1213790)



Catalog Number: 176896

Product name

Osocimab (Synonyms: BAY1213790)

Specificity

Human F11 / Factor XI

Antibody description

Anti-F11 / Factor XI Antibody (Osocimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

1×PBS, pH 6.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG

Applications


ELISA

Validations



SDS-PAGE

Anti-F11 / Factor XI Antibody (osocimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



SEC-HPLC

The purity of Anti-F11 / Factor XI Antibody (osocimab) is more than 95%, determined by SEC-HPLC.