

# Trevogrumab (Synonyms: HUMAN IGG4 ANTI-MYOSTATIN MAB)



Catalog Number: 176912

## Product name

IgG4

Trevogrumab (Synonyms: HUMAN IGG4 ANTI-MYOSTATIN MAB)


## Applications

ELISA

## Specificity

Human GDF8 / Myostatin

## Validations



## Antibody description


Anti-GDF8 / Myostatin Antibody (Trevogrumab)

## Preparation

Recombinant expression and purified from CHO cells.

## SDS-PAGE

Anti-GDF8 / Myostatin Antibody (trevogrumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## SEC-HPLC

The purity of Anti-GDF8 / Myostatin Antibody (trevogrumab) is 96.57%, determined by SEC-HPLC.

## Clonality

Monoclonal

## Ig Type