

Amlitelimab reference antibody



Catalog Number: 176927

Product name

Amlitelimab reference antibody

Specificity

Human TNFSF4 / OX40L / CD252

Antibody description

Anti-TNFSF4 / OX40L / CD252 Antibody (Amlitelimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG4SP

Applications


ELISA

Validations



SDS-PAGE

Anti-TNFSF4 / OX40L / CD252 Antibody (amlitelimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



SEC-HPLC

The purity of Anti-TNFSF4 / OX40L / CD252 Antibody (amlitelimab) is more than 95%, determined by SEC-HPLC.