

Lumiliximab (Synonyms: PRIMATIZED ANTI-CD23; GOMILIXIMAB [USAN])



Catalog Number: 176932

Product name

IgG1

Lumiliximab (Synonyms: PRIMATIZED ANTI-CD23; GOMILIXIMAB [USAN])


Applications

ELISA

Specificity

Human FcεR2 / CD23

Validations



Antibody description


Anti-FcεR2 / CD23 Antibody (Lumiliximab)

Preparation

Recombinant expression and purified from CHO cells.

SDS-PAGE

Anti-FcεR2 / CD23 Antibody (lumiliximab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

SEC-HPLC

The purity of Anti-FcεR2 / CD23 Antibody (lumiliximab) is 100%, determined by SEC-HPLC.

Clonality

Monoclonal

Ig Type