

Tesnatilimab reference antibody



Catalog Number: 176940

Product name

Tesnatilimab reference antibody

Specificity

Human NKG2D / CD314

Antibody description

Anti-NKG2D / CD314 Antibody (Tesnatilimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG4

Applications


ELISA

Validations




SDS-PAGE

Anti-NKG2D / CD314 Antibody (tesnatilimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-NKG2D / CD314 Antibody (tesnatilimab) is 100%, determined by SEC-HPLC.



Bioactivity: ELISA

Immobilized human NKG2D His at 2 µg/mL can bind Anti-NKG2D / CD314 Antibody (tesnatilimab), EC50=0.0209 µg/mL.