

Camrelizumab (Synonyms: Carilizumab; Carrelizumab)



Catalog Number: 176941

Product name

Camrelizumab (Synonyms: Carilizumab; Carrelizumab)

Specificity

Human PDCD1 / PD-1 / CD279

Antibody description

Anti-PDCD1 / PD-1 / CD279 Antibody (Camrelizumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG4SP

Applications


ELISA

Validations




SDS-PAGE

Anti-PDCD1 / PD-1 / CD279 Antibody (camrelizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-PDCD1 / PD-1 / CD279 Antibody (camrelizumab) is 100%, determined by SEC-



HPLC

Bioactivity: ELISA

Immobilized human PD-1 His at 2 ug/mL can bind Anti-PDCD1 / PD-1 / CD279 Antibody (camrelizumab), EC50=0.008247 ug/mL.