

# Vibostolimab reference antibody



Catalog Number: 176950

## Product name

Vibostolimab reference antibody

## Specificity

Human TIGIT

## Antibody description

Anti-TIGIT Antibody (Vibostolimab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-TIGIT Antibody (vibostolimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-TIGIT Antibody (vibostolimab) is 100%, determined by SEC-HPLC.



Bioactivity: ELISA

Immobilized human TIGIT His at 2 µg/mL can bind Anti-TIGIT Antibody (vibostolimab), EC50=0.01861 µg/mL.