

# Quetmolimab reference antibody



Catalog Number: 176952

## Product name

Quetmolimab reference antibody

## Specificity

Human CX3CL1 / Fractalkine

## Antibody description

Anti-CX3CL1 / Fractalkine Antibody  
(Quetmolimab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG2SA

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-CX3CL1 / Fractalkine Antibody (quetmolimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-CX3CL1 / Fractalkine Antibody (quetmolimab) is 98.99%, determined by SEC-HPLC.