

Begelomab (Synonyms: BEGEDINA)



Catalog Number: 176953

Product name

Begelomab (Synonyms: BEGEDINA)

Specificity

Human DPP4 / CD26

Antibody description

Anti-DPP4 / CD26 Antibody (Begelomab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG2SA

Applications


ELISA

Validations



SDS-PAGE

Anti-DPP4 / CD26 Antibody (begelomab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-DPP4 / CD26 Antibody (begelomab) is 98.65%, determined by SEC-HPLC.