

# Etrolizumab (Synonyms: ANTI-.BETA.7)



Catalog Number: 176958

## Product name

Etrolizumab (Synonyms: ANTI-.BETA.7)

## Specificity

Human Integrin a4b7 (ITGA4 & ITGB7)

## Antibody description

Anti-Integrin a4b7 (ITGA4 & ITGB7) Antibody (Etrolizumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-Integrin a4b7 (ITGA4 & ITGB7) Antibody (etrolizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-Integrin a4b7 (ITGA4 & ITGB7) Antibody (etrolizumab) is 100%, determined by SEC-HPLC.