

# Revdofilimab reference antibody



Catalog Number: 176971

## Product name

Revdofilimab reference antibody

## Specificity

Human TNFRSF4 / OX40 / CD134

## Antibody description

Anti-TNFRSF4 / OX40 / CD134 Antibody  
(Revdofilimab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-TNFRSF4 / OX40 / CD134 Antibody (revdofilimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-TNFRSF4 / OX40 / CD134 Antibody (revdofilimab) is more than 95%, determined by SEC-HPLC.