

# Bersanlimab reference antibody



Catalog Number: 176981

## Product name

Bersanlimab reference antibody

## Specificity

Human ICAM1 / CD54

## Antibody description

Anti-ICAM1 / CD54 Antibody (Bersanlimab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations




## SEC-HPLC

The purity of Anti-ICAM1 / CD54 Antibody (bersanlimab) is 100%, determined by SEC-HPLC.



## ELISA

Immobilized human ICAM-1 / CD54 Protein, His Tag Protein at 2 ug/mL can bind anti-ICAM1 / CD54 Antibody (bersanlimab)  $EC_{50}=0.002577$



ug/mL.

## SDS-PAGE

Anti-ICAM1 / CD54 Antibody (bersanlimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95.0%.