

Ziltivekimab (Synonyms: MEDI 5117 [WHO-DD])



Catalog Number: 176983

Product name

Ziltivekimab (Synonyms: MEDI 5117 [WHO-DD])

Specificity

Human IL-6 / IFN γ 2

Antibody description

Anti-IL-6 / IFN γ 2 Antibody (Ziltivekimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

1 \times PBS, pH 6.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG

Applications


ELISA

Validations



SDS-PAGE

Anti-IL-6 / IFN γ 2 Antibody (ziltivekimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



SEC-HPLC

The purity of Anti-IL-6 / IFN γ 2 Antibody (ziltivekimab) is more than 95%, determined by SEC-HPLC.