

Belantamab reference antibody



Catalog Number: 176989

Product name

Belantamab reference antibody

Specificity

Human TNFRSF17 / BCMA / CD269

Antibody description

Anti-TNFRSF17 / BCMA / CD269 Antibody
(Belantamab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations




SDS-PAGE

Anti-TNFRSF17 / BCMA / CD269 Antibody (belantamab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-TNFRSF17 / BCMA / CD269 Antibody (belantamab) is 100%, determined by



SEC-HPLC.

Bioactivity: FACS

U266 cells were stained with Anti-TNFRSF17 / BCMA / CD269 Antibody (belantamab) and negative control protein respectively, washed and then followed by PE and analyzed with FACS, EC50=1.3980 µg/mL.