

# Ascrinvacumab reference antibody



Catalog Number: 176993

## Product name

Ascrinvacumab reference antibody

## Specificity

Human ACVRL1 / ALK-1

## Antibody description

Anti-ACVRL1 / ALK-1 Antibody (Ascrinvacumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

1×PBS, pH 6.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-ACVRL1 / ALK-1 Antibody (ascrinvacumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



## SEC-HPLC

The purity of Anti-ACVRL1 / ALK-1 Antibody (ascrinvacumab) is more than 95%, determined by SEC-HPLC.