

Grisnilimab reference antibody



Catalog Number: 177024

Product name

Grisnilimab reference antibody

Specificity

Human CD7

Antibody description

Anti-CD7 Antibody (Grisnilimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG2SA

Applications


ELISA

Validations



SDS-PAGE

Anti-CD7 Antibody (grisnilimAb) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



SEC-HPLC

The purity of Anti-CD7 Antibody (grisnilimAb) is 94.62%, determined by SEC-HPLC.