

Sirtratumab reference antibody



Catalog Number: 177035

Product name

Sirtratumab reference antibody

Specificity

Human SLITRK6

Antibody description

Anti-SLITRK6 Antibody (Sirtratumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG2SA

Applications


ELISA

Validations




SDS-PAGE

Anti-SLITRK6 Antibody (sirtratumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-SLITRK6 Antibody (sirtratumab) is 100%, determined by SEC-HPLC.



Bioactivity: ELISA

Immobilized human SLITRK6 His at 2 µg/mL can bind Anti-SLITRK6 Antibody (sirtratumab), EC50=0.04938 µg/mL.