

Ordesekimab (Synonyms: AMG-714; HUMAX-IL 15)



Catalog Number: 177071

Product name

IgG1

Ordesekimab (Synonyms: AMG-714; HUMAX-IL 15)


Applications

ELISA

Specificity

Human IL-15

Validations



Antibody description

Anti-IL-15 Antibody (Ordesekimab)

Preparation

Recombinant expression and purified from CHO cells.

SDS-PAGE


Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Anti-IL-15 Antibody (ordesekimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.



Clonality

Monoclonal

SEC-HPLC

Ig Type

The purity of Anti-IL-15 Antibody (ordesekimab) is 95.77%, determined by SEC-HPLC.