

Tadocizumab reference antibody



Catalog Number: 177093

Product name

Tadocizumab reference antibody

Specificity

Human Integrin $\alpha 2\beta 3$ (ITGA2 & ITGB3)

Antibody description

Anti-Integrin $\alpha 2\beta 3$ (ITGA2 & ITGB3) Antibody (Tadocizumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Glycine, 0.2 NaCl, 50 mM Tris-HCl, pH 6.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations



SDS-PAGE

Anti-Integrin $\alpha 2\beta 3$ (ITGA2 & ITGB3) Antibody (tadocizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



SEC-HPLC

The purity of Anti-Integrin $\alpha 2\beta 3$ (ITGA2 & ITGB3) Antibody (tadocizumab) is more than 95%, determined by SEC-HPLC.