

# Mogamulizumab reference antibody



Catalog Number: 177110

## Product name

Mogamulizumab reference antibody

## Specificity

Human CCR4 / CD194

## Antibody description

Anti-CCR4 / CD194 Antibody (Mogamulizumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

1×PBS, pH 6.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-CCR4 / CD194 Antibody (mogamulizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



## SEC-HPLC

The purity of Anti-CCR4 / CD194 Antibody (mogamulizumab) is more than 95%, determined by SEC-HPLC.