

# Tinurilimab reference antibody



Catalog Number: 177114

## Product name

Tinurilimab reference antibody

## Specificity

Human CEACAM6 / CD66c

## Antibody description

Anti-CEACAM6 / CD66c Antibody (Tinurilimab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG2SA

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-CEACAM6 / CD66c Antibody (tinurilimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-CEACAM6 / CD66c Antibody (tinurilimab) is 95.93%, determined by SEC-HPLC.



## Bioactivity: ELISA

Immobilized human CEACAM6 His at 2 µg/mL can bind Anti-CEACAM6 / CD66c Antibody (tinurilimab), EC<sub>50</sub>=0.005084 µg/mL.