

# Enoblituzumab reference antibody



Catalog Number: 177128

## Product name

Enoblituzumab reference antibody

## Specificity

Human B7-H3 / CD276

## Antibody description

Anti-B7-H3 / CD276 Antibody (Enoblituzumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-B7-H3 / CD276 Antibody (enoblituzumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-B7-H3 / CD276 Antibody (enoblituzumab) is 98.61%, determined by SEC-



## HPLC.

Bioactivity: ELISA

Immobilized human B7H3 His at 2 ug/mL can bind Anti-B7-H3 / CD276 Antibody (enoblituzumab), EC50=0.003300 ug/mL.