

Gatipotuzumab reference antibody



Catalog Number: 177141

Product name

Gatipotuzumab reference antibody

Monoclonal

Ig Type

IgG1

Specificity

Human MUC1


Applications

ELISA

Antibody description

Anti-MUC1 Antibody (Gatipotuzumab)

Validations



Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

SDS-PAGE

Anti-MUC1 Antibody (gatipotuzumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality