

Gatipotuzumab reference antibody



Catalog Number: 177141

Product name

Monoclonal

Gatipotuzumab reference antibody

Ig Type

Specificity

IgG1

Human MUC1

Applications

Antibody description

ELISA

Anti-MUC1 Antibody (Gatipotuzumab)

Validations

Preparation

Recombinant expression and purified from CHO cells.


Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality



SDS-PAGE

Anti-MUC1 Antibody (gatipotuzumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.